Varicex® S

Bi-elastische Zinkleimbinde

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2. Usage :
 - Varicex® can be applied using a variety of techniques. The bandaging technique is chosen in consultation with the treating physician.
 - To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
 - Apply the bandage so that there are no gaps.
 - The bandage is unrolled with the rolled-up portion (roll) on top. Keep the roll close to the leg / arm.
 - The bandage must be replaced when it becomes loose.

Due to its bi-elasticity, Varicex® S can be moulded to the body without making multiple clips or cuts (Fig. 1-6). The bandage conforms to the leg without creases.
Varicex® F – if necessary to prevent gaps in the dressing, cut off the bandage and then start wrapping it again (Fig. 10 and 11).

After it has been applied, the zinc paste bandage should be protected by a tg® Tubular Bandage or a conforming bandage (Fig. 9-11).

Disposal:
The product must be disposed of as hazardous waste, since it contains zinc oxide, which is classified as chronically hazardous to the aquatic environment.
The user must observe the respective laws, standards and guidelines on disposal of waste from the medical sector.
In Europe, non-contaminated packaging waste can be assigned a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

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Dispositif médical

en Instructions for use

Product description and performance characteristics:
Varicex® S is a short-stretch, bi-elastic zinc paste bandage. Suitable for an Unna boot.
Varicex® F is a rigid zinc past bandage. Suitable for an Unna boot.

Product composition:
Varicex® S: backing material: 80% viscose, 20% polyamide - with zinc paste coating
Varicex® F: backing material: 40% cotton, 60% viscose - with zinc paste coating
Zinc paste: zinc oxide, gelatine, calcium chloride, glycerol, water, preservatives (propyl 4-hydroxybenzoate, methyl 4-hydroxybenzoate) and fragrance.

Intended purpose:
For compression, support and relief, for immobilisation of limbs.

Indications:
Phlebological and lymphological indications: varicosis, deep venous insufficiency, thrombophlebitis, phlebotrombosis, post-thrombotic syndrome, Widmer's stage I to III b chronic venous insufficiency, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery.
Contusions, sprains, tendon injuries and joint dislocations.

Contraindications:
Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases, skin irritations and/or allergies may occur.

- Warnings and precautions:**
- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
 - Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk. If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression.
 - The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
 - For use with ABPIs over 0.8. When the ABPI is between 0.5 and 0.8, apply under medical supervision only.
 - Known allergy and/or hypersensitivity to any of the product components.

- Usage information:**
1. Preparation:
 - Before applying the bandage, clean / debride any wound that is present and cover with a wound dressing appropriate for the wound status.
 - If necessary, apply padding to areas at risk of developing pressure ulcers.
 - It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
 - The width of the bandage should correspond to the diameter of the limb to be banded.
 2. Usage:
 - Varicex® can be applied using a variety of techniques. The bandaging technique is chosen in consultation with the treating physician.
 - To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
 - Apply the bandage so that there are no gaps.
 - The bandage is unrolled with the rolled-up portion (roll) on top. Keep the roll close to the leg / arm.
 - The bandage must be replaced when it becomes loose.

Due to its bi-elasticity, Varicex® S can be moulded to the body without making multiple clips or cuts (Fig. 1-6). The bandage conforms to the leg without creases.
Varicex® F – if necessary to prevent gaps in the dressing, cut off the bandage and then start wrapping it again (Fig. 10 and 11).

After it has been applied, the zinc paste bandage should be protected by a tg® Tubular Bandage or a conforming bandage (Fig. 9-11).

Disposal:
The product must be disposed of as hazardous waste, since it contains zinc oxide, which is classified as chronically hazardous to the aquatic environment.
The user must observe the respective laws, standards and guidelines on disposal of waste from the medical sector.
In Europe, non-contaminated packaging waste can be assigned a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

Note:
Contains zinc oxide, avoid release into the environment (poisonous for water organisms, chronically hazardous to the aquatic environment), please take note of the special information in the "Disposal" section.

General instructions:
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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