

PRODUCT SPECIFICATION

Product name:

sempercare® green

Spec code:

NCF-035GR-N-3CZ

Date of issue:

14-Nov-2017

Type Labelling Shape Material Colour Inner Outer Cuff/surface Storage life

single use examination and disposable protective glove, non sterile all information on dispenser box ambidextrous - straight fingers Nitrile Butadiene Rubber (NBR) - [accelerator-free and latex-free] green

polymer coated / powder free no treatment rolled cuff/fingertip textured 3 years

Available sizes

XS (5-6)

S (6-7) M (7-8)

L (8-9)

XL (9-10)

Min. thickness [in mm]	Finger	Palm	Cuff
single measured	0.07	0.06	0.05
double measured	0.14	0.12	0.10

	Before aging	After aging	Testing method
Elongation	min. 500 %	min. 400 %	ASTM D6319
Tensile strength	min. 14 MPa	min. 14 MPa	ASTM D6319

AQL (for pinholes)	≤ 1.5	EN 374-2 EN 455-1
Residual powder / Powder content	≤ 2 mg/glove	EN 455-3
Force at break	median ≥ 6.0 N	EN 455-2
Length	median 240 mm	EN 455-2

Storage instructions

Store according to ISO 2230:2002. Ensure that the area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

Page 1/2

sempercare® green

NCF-035GR-N-3CZ_Version A_November 2017_EN

Replaces all previous versions and valid for 2 years after date of issue

Standards, guidelines & quality certificates

Quality certification	ISO 9001 ISO 13485 ISO 14001 OSHAS 18001
Conformity to directives	Personal Protective Equipment 89/686/EEC: Cat III for complex risks
	Medical Device Directive 93/42/EEC: Class I
Confomity to standards	EN 420 EN 455 1-4 EN ISO 374-1, EN 374-2, ISO 16523, EN 374-4
	ASTM D6319 ASTM F1671
Conformity to food regulation	1935/2004 (EC)
labelling according to	ISO 15223 EN 420 EN ISO 374

Warning and reporting system

Potentially allergenic ingredients

not known

For further information, the actual version of the ingredients list is given on request.

Medical device vigilance and reporting system

According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. E-Mail: sempermed.complaints@semperitgroup.com or Tel.: +43 2630 310 ext. -510; Fax: -549

Authorized signatory

C. Eichler

Product management

L. Millonig

Regulatory affairs

J. Glantschnig

Page 2/2

sempercare® green

NCF-035GR-N-3CZ_Version A_November 2017_EN

Replaces all previous versions and valid for 2 years after date of issue