




|   |  |           |   |          |   |                    |
|---|--|-----------|---|----------|---|--------------------|
|   | NATURAL RUBBER LATEX (NR)  |           |   |          |   |                    |
| BRAND NAME  | sempermed® supreme   |           | sempermed® supreme plus   |          | sempermed® supreme green  |                    |
| <div><div><div>quick don™</div></div><div><div>quick don™</div></div><div><div>quick don™</div></div></div> |  |           |   |          |   |                    |
| PRODUCT BENEFITS  | <b>A PREMIUM ALL-ROUNDER YOU CAN TRUST</b> <ul style="list-style-type: none"><li>■ <b>Extra-supple:</b><br/>great comfort &amp; tactile sensitivity</li><li>■ <b>Fast &amp; easy to don:</b><br/>Quick Don™ inner coating</li><li>■ <b>Anatomical shape:</b><br/>reduces strain for wearer</li><li>■ <b>Skin-friendly:</b><br/>washed to minimise skin irritants</li></ul> |           | <b>OPTIMAL GRIP IN EVERY SITUATION</b> <ul style="list-style-type: none"><li>■ <b>Micro-rough surface:</b><br/>for dry &amp; damp applications</li><li>■ <b>Fast &amp; easy to don:</b><br/>Quick Don™ inner coating</li><li>■ <b>Anatomical shape:</b><br/>reduces strain for wearer</li><li>■ <b>Skin-friendly:</b><br/>washed to minimise skin irritants</li></ul> |          | <b>DOUBLE PROTECTION THAT YOU CAN TRUST</b> <ul style="list-style-type: none"><li>■ <b>Clear detection of perforations:</b><br/>when double-donning with a lighter coloured outer glove</li><li>■ <b>Tactile sensitivity</b><br/>thanks to optimised wall thickness</li><li>■ <b>Ideal inner glove:</b><br/>designed for double-donning</li><li>■ <b>As underglove</b> for sempermed® supreme or supreme plus</li></ul> |                    |
| TYPE  | Sterile surgical glove for single use, powder free, with synthetic lining  |           |   |          |   |                    |
| MATERIAL  | Natural rubber latex (NR)  |           | Natural rubber latex (NR)   |          | Natural rubber latex (NR)   |                    |
| COLOUR  | Natural white  |           | Natural white   |          | Green   |                    |
| GLOVE SHAPE   | Fully anatomical with curved fingers and rolled rim  |           |   |          |   |                    |
| SIZE/OVERALL LENGTH<br>as per EN 455-2  |  |           | 5 ½, 6 and 6 ½    median 270 mm<br>7, 7 ½ and 8      median 280 mm<br>8 ½ and 9          median 285 mm  |          |   |                    |
| WALL THICKNESS measured<br>in single layer in the palm area   | 0.19–0.23 mm   |           | 0.19–0.23 mm  |          | 0.17–0.20 mm  |                    |
| BARRIER PERFORMANCE   | AQL 0.65   |           |   |          |   |                    |
| FORCE AT BREAK<br>as per EN 455-2   | ≥ 12 N, exceeding EN 455-2 (≥ 12 N instead of 9 N)   |           |   |          |   |                    |
| DURABILITY<br>in original packaging if stored<br>as per DIN 7716, ISO 2230  | 3 years  |           |   |          |   |                    |
| STERILISATION   | Radiation <div>STERILE</div> <div>R</div>  |           |   |          |   |                    |
| PACKAGING   | Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack.<br>In dispenser carton with sterilisation indicator: 50 pairs.<br>In transport carton with sterilisation indicator: 300 pairs.   |           |   |          |   |                    |
| PURPOSE according to MP<br>Regulation 93/42/EWG / according<br>to PPE Regulation (EU) 2016/425  | Medical device class IIa<br>Single-use protective glove – PPE category III<br>(protection against chemical substances for limited time) <sup>2)</sup>  |           |   |          |   |                    |
| ARTICLE<br>NUMBERS  | EU   |           | INT   |          | EU  |                    |
|   | Size 5 ½   | 822751521 | 822751525   | Size 5 ½ | 822851521   | Size 5 ½ 823756521 |
|   | Size 6   | 822751601 | 822751605   | Size 6   | 822851601   | Size 6 823756601   |
|   | Size 6 ½   | 822751621 | 822751625   | Size 6 ½ | 822851621   | Size 6 ½ 823756621 |
|   | Size 7   | 822751701 | 822751705   | Size 7   | 822851701   | Size 7 823756701   |
|   | Size 7 ½   | 822751721 | 822751725   | Size 7 ½ | 822851721   | Size 7 ½ 823756721 |
|   | Size 8   | 822751801 | 822751805   | Size 8   | 822851801   | Size 8 823756801   |
|   | Size 8 ½   | 822751821 | 822751825   | Size 8 ½ | 822851821   | Size 8 ½ 823756821 |
|   | Size 9   | 822751901 | 822751905   | Size 9   | 822851901   | Size 9 823756901   |

<sup>1)</sup> Crepy et.al. (2018), Accelerator-free gloves as alternatives in cases of glove allergy in healthcare workers, in Contact Dermatitis. 2018 Jan;78(1):28-32. doi: 10.1111/cod.12860. Epub 2017 Jul 27. <https://www.ncbi.nlm.nih.gov/pubmed/28748553>

<sup>2)</sup> Please note the Sempermed Chemical Resistance list, available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>