

BRAND NAME

**sempermed<sup>®</sup>**  
**syntegra UV**

**sempermed<sup>®</sup>**  
**syntegra IR**

**sempermed<sup>®</sup>**  
**syntegra green**



**PRODUCT BENEFITS**

**AWARD-WINNING PROTECTION FOR YOUR SKIN**

- **Polyisoprene:**  
excellent comfort & protection
- **Low-allergenic<sup>1)</sup>:**  
accelerator-free & latex free
- **Great comfort:**  
extraordinary latex-like fit & feel
- **Fast & easy to don:**  
Quick Don™ inner coating

**ADVANCED SKIN PROTECTION & PERFORMANCE**

- **Reduced allergy risk:**  
innovative accelerators, latex free
- **Fast & easy to don:**  
Quick Don™ inner coating
- **Great comfort:**  
extraordinary latex-like fit & feel
- **Anatomical shape:**  
reduces strain for wearer

**THE INNER GLOVE FOR HIGH-RISK PROCEDURES**

- **Ideal inner glove:**  
designed for double-donning
- **Skin-friendly:**  
innovative accelerators, latex free
- **Great comfort:**  
extraordinary latex-like fit & feel
- **Wall thickness**  
optimal for high-risk surgeries

**TYPE** Sterile surgical glove for single use, powder free, with synthetic lining

**MATERIAL** Synthetic polyisoprene Synthetic polyisoprene Synthetic polyisoprene

**COLOUR** White Cream Green

**GLOVE SHAPE** Fully anatomical with curved fingers and rolled rim

**SIZE/OVERALL LENGTH**  
as per EN 455-2

5 ½, 6 and 6 ½	median 270 mm
7, 7 ½ and 8	median 280 mm
8 ½ and 9	median 285 mm

**WALL THICKNESS** measured in single layer in the palm area

0.19–0.25 mm	0.19–0.24 mm	0.19–0.24 mm
--------------	--------------	--------------

**BARRIER PERFORMANCE** AQL 0.65

**FORCE AT BREAK** as per EN 455-2 ≥ 9 N

**DURABILITY** in original packaging if stored as per DIN 7716, ISO 2230 3 years

**STERILISATION** Radiation STERILE R

**PACKAGING**

Left and right glove turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 40 pairs. In transport carton with sterilisation indicator: 240 pairs.	Left and right glove turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs. In transport carton with sterilisation indicator: 300 pairs.
---	---

**PURPOSE** according to MP Regulation 93/42/EWG / according to PPE Regulation (EU) 2016/425

Medical device class IIa  
Single-use protective glove – PPE category III  
(protection against chemical substances for limited time)<sup>2)</sup>



ARTICLE NUMBERS	EU		EU		EU	
	Size	Article No.	Size	Article No.	Size	Article No.
	Size 5 ½	827058521	Size 5 ½	827056521	Size 5 ½	827053521
	Size 6	827058601	Size 6	827056601	Size 6	827053601
	Size 6 ½	827058621	Size 6 ½	827056621	Size 6 ½	827053621
	Size 7	827058701	Size 7	827056701	Size 7	827053701
	Size 7 ½	827058721	Size 7 ½	827056721	Size 7 ½	827053721
	Size 8	827058801	Size 8	827056801	Size 8	827053801
	Size 8 ½	827058821	Size 8 ½	827056821	Size 8 ½	827053821
	Size 9	827058901	Size 9	827056901	Size 9	827053901

<sup>1)</sup> Crepy et.al. (2018), Accelerator-free gloves as alternatives in cases of glove allergy in healthcare workers, in Contact Dermatitis. 2018 Jan;78(1):28-32. doi: 10.1111/cod.12860. Epub 2017 Jul 27. <https://www.ncbi.nlm.nih.gov/pubmed/28748553>

<sup>2)</sup> Please note the Sempermed Chemical Resistance list, available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>